## BIOLOGY FORM 4 CONFIDENTIAL

About 3ml of solution L- enzyme diastase per student About 3ml of 0.05% NaCl solution with a dropper

About 3ml of 1.4% NaCl solution with a dropper

About 10ml of starchsolution per student

3 test tubes

Measuring cylinder (5 ml or 10 ml)

White tile

Iodine solution with a dropper

Benedict's solution with a dropper

Test tube holder

Warm water bath maintained between 35-38°C

Thermometer

Means of labeling for the three test tubes

1 dropper

Test tube rack

Hibiscus leaf labelled P

Tradescantia leaf labelled R

Bougainvillea twig with a flower labelled T

Maize/ Grass leaf labelled X

Bidens pilosa leaf labelled Y